

## Product Datasheet

### Secretory Immunoglobulin A, Human Colostrum 527-55

<b>Product Name</b>	Secretory Immunoglobulin A, Human Colostrum
<b>Catalog Number</b>	527-55
<b>Alternative Names</b>	slgA
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	Reported in mg/vial ( $A_{278} E^{0.1\%}=1.26$ )
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Human Colostrum
<b>Cross Reactivity</b>	Non-reactive with antisera to human albumin, IgD, IgE, IgG, IgM, and lactoferrin (IEP)
<b>Molecular Weight</b>	~385,000
<b>Shelf Life</b>	3 years