

## Product Datasheet

### Anti-Human CD27 Monoclonal Antibody, Unconjugated Clone LT27 853.343.020

**Product Name** Anti-Human CD27 Monoclonal Antibody, Unconjugated Clone LT27

**Catalog Number** 853.343.020

**Description** Mouse monoclonal

**Species/Host** Mouse

**Conjugation** Purified

**Tested Applications** FC

**Immunogen** Human peripheral blood lymphocytes

**Brand** Diaclone

**Form/Appearance** Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4

**Storage** Store at +2-8°C. DO NOT FREEZE.

**Note** We recommend that each laboratory determine an optimum working titre for use in its particular application

**Isotype** IgG2a

**Clonality** Monoclonal

**Specificity** Reacts with CD27 (T14), a 50-55 kDa type I transmembrane glycoprotein (member of the TNF-receptor superfamily)