

Product Datasheet

Human Heart Troponin Complex (I-T-C) 550-10

Product Name	Human Heart Troponin Complex (I-T-C)
Catalog Number	550-10
Alternative Names	(I-T-C), Blend of Troponin I, Troponin C and Troponin T
Brand	Lee Biosolutions
Form/Appearance	Liquid Frozen, brown liquid
Concentration	Report mg/mL (Lowry)
Storage	-60 °C or below
Purity	Partially Purified
CAS Number	77108-40-8
Source	Human Heart
Shelf Life	5 years
Biological Activity	Troponin I: $\geq 25,000$ ng/mL (Immunoassay)