

Product Datasheet

Troponin I (Tn-I), Human Heart 550-11

Product Name	Troponin I (Tn-I), Human Heart
Catalog Number	550-11
Alternative Names	Tn-I, hTn-I, cTnI, cardiac Tn-I, TnI
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.5 mg/mL (A_{280} , $E^{0.1\%}=0.42$)
Storage	-20 °C
Purity	$\geq 98\%$ (SDS-PAGE)
CAS Number	77108-40-8
Source	Human Heart
Associated Products	Anti-h cTnI 9701/RC970, Anti-h cTnI 9703, Anti-h cTnI 9705, Anti-h cTnI 9707/RC9707, Anti-h cTnI RC9750
Molecular Weight	~24,000
Shelf Life	2 years