

Product Datasheet

Troponin T (Tn-T), Human Heart 550-21

Product Name	Troponin T (Tn-T), Human Heart
Catalog Number	550-21
Alternative Names	Tn-T, hTn-T, cardiac Tn-T, TnT, cTnT, cTn-T
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.5 mg/mL (A ₂₈₀)
Storage	-20 °C (up to three months), -70 °C (more than three months)
Purity	≥ 98% (SDS-PAGE)
CAS Number	77108-40-8
Source	Human Heart
Molecular Weight	~35,000
Shelf Life	2 years