

Product Datasheet

Breast Milk, Pooled Human Donors - ≥4 Weeks 991-01-P

Product Name	Breast Milk, Pooled Human Donors - ≥4 Weeks
Catalog Number	991-01-P
Alternative Names	Lactate, Lactation
Brand	Lee Biosolutions
Form/Appearance	Liquid White to cream colored liquid
Storage	-20 °C
Source	Pooled Human Donors
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None