

## Product Datasheet

### Breast Milk, Pooled Human Donors - ≥4 Weeks 991-01-P

<b>Product Name</b>	Breast Milk, Pooled Human Donors - ≥4 Weeks
<b>Catalog Number</b>	991-01-P
<b>Alternative Names</b>	Lactate, Lactation
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid White to cream colored liquid
<b>Storage</b>	-20 °C
<b>Source</b>	Pooled Human Donors
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
<b>Preservatives</b>	None