

## Product Datasheet

### Pooled Human Urine - Normal - Fresh, Never Frozen 991-03-P-FNF

<b>Product Name</b>	Pooled Human Urine - Normal - Fresh, Never Frozen
<b>Catalog Number</b>	991-03-P-FNF
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Light to dark yellow, clear liquid
<b>Storage</b>	2-8° C
<b>Purity</b>	Raw Urine
<b>Source</b>	Pooled Human Donors
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
<b>Preservatives</b>	None (Available by request)