

## Product Datasheet

### Pooled Human Urine - pH Treated 991-03-P-PH

**Product Name** Pooled Human Urine - pH Treated

**Catalog Number** 991-03-P-PH

**Brand** Lee Biosolutions

**Form/Appearance** Liquid

**Storage** -20 °C

**Purity** pH Treated Urine

**Source** Pooled Human Donors

**Shelf Life** Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.

**Preservatives** None