

Product Datasheet

Anti-Human CD35 Monoclonal Antibody, Unconjugated Clone UJ11 853.403.020

Product Name Anti-Human CD35 Monoclonal Antibody, Unconjugated Clone UJ11

Catalog Number 853.403.020

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications FC

Brand Diaclone

Form/Appearance PBS with 15 mM sodium azide (Ph7.4)

Storage Store at 4°C. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.

Note We recommend that each laboratory determine an optimum working titre for use in its particular application

Isotype IgG1,k

Clonality Monoclonal

Specificity The CD35 antigen, CR1, is a complement component receptor with specificity for C3b, C4b and iC3b