

Product Datasheet

Anti-Human CD57 Monoclonal Antibody, Unconjugated Clone NC1 853.463.020

Product Name Anti-Human CD57 Monoclonal Antibody, Unconjugated Clone NC1

Catalog Number 853.463.020

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications IHC, FC

Immunogen Quail ciliary ganglia cells

Brand Diaclone

Form/Appearance 20 mM TRIS and 150 mM NaCl with 1% BSA, pH 8.0. Preservative: Mergiolate (0.001%)

Storage Store at +2-8°C. For longer storage freeze aliquots.

Note We recommend that each laboratory determine an optimum working titre for use in its particular application

Isotype IgM

Clonality Monoclonal

Source Quail ciliary ganglia cells

Specificity Directed against the CD57-antigen, which is expressed on NK-cells and subset T lymphocytes

Hybridoma Myeloma SP2/O Ag 1.4 x Balb/c spleen cells