

## Product Datasheet

### Human anti-Mouse Antibody (HAMA) Serum 991-24-S-HAMA

<b>Product Name</b>	Human anti-Mouse Antibody (HAMA) Serum
<b>Catalog Number</b>	991-24-S-HAMA
<b>Alternative Names</b>	HAMA Serum
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20 °C
<b>Source</b>	Single Human Donor
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.