

## Product Datasheet

### High NT-proBNP Serum, Human Donor 991-24-S-NTPBNP

<b>Product Name</b>	High NT-proBNP Serum, Human Donor
<b>Catalog Number</b>	991-24-S-NTPBNP
<b>Tested Applications</b>	Control
<b>Alternative Names</b>	NT-Pro Brain Natriuretic Peptide, Brain Natriuretic Peptide; proBNP; Natriuretic Peptides, B-type Natriuretic Peptide; N-terminal pro b-type Natriuretic Peptide
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20 °C
<b>Source</b>	Single Human Donor
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.