

## Product Datasheet

### Perianal Swabs - Human Female Donor 991-34-F

|                          |  |
|--------------------------|--|
| <b>Product Name</b>      | Perianal Swabs - Human Female Donor  |
| <b>Catalog Number</b>    | 991-34-F   |
| <b>Alternative Names</b> | Perineum swabs   |
| <b>Brand</b>             | Lee Biosolutions   |
| <b>Form/Appearance</b>   | Swab   |
| <b>Storage</b>           | -20 °C   |
| <b>Source</b>            | Single Human Donor   |
| <b>Shelf Life</b>        | Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time. |