

Product Datasheet

Immunoglobulin A (IgA), anti-Human Polyclonal Antibody GA-80A

Product Name Immunoglobulin A (IgA), anti-Human Polyclonal Antibody

Catalog Number GA-80A

Species/Host Goat

Immunogen High purity native human plasma Immunoglobulin A

Alternative Names anti-IgA

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration 1 - 2 mg/mL ($E^{1\text{mg/mL}}_{280\text{nm}} = 1.4$)

Storage 2-8° C

Isotype IgG

Purity Affinity Purified

Source Polyclonal Antibody

Shelf Life 2 years