

## Product Datasheet

### C4 Polyclonal Goat Anti-Human Affinity Purified Antibody GC4-80A

**Product Name** C4 Polyclonal Goat Anti-Human Affinity Purified Antibody

**Catalog Number** GC4-80A

**Species/Host** Goat

**Immunogen** High purity native human plasma Complement 4 (C4)

**Brand** Lee Biosolutions

**Form/Appearance** Liquid Clear solution

**Storage** 2-8° C

**Isotype** IgG

**Purity** Affinity Purified

**Source** Polyclonal Antibody

**Shelf Life** 2 years