
Product Datasheet

C4 Polyclonal Goat Anti-Human Affinity Purified Antibody GC4-80A

Product Name C4 Polyclonal Goat Anti-Human Affinity Purified Antibody

Catalog Number GC4-80A

Species/Host Goat

Immunogen High purity native human plasma Complement 4 (C4)

Brand Lee Biosolutions

Form/Appearance Liquid Clear solution

Storage 2-8° C

Isotype IgG

Purity Affinity Purified

Source Polyclonal Antibody

Shelf Life 2 years