

## Product Datasheet

### C-Reactive Protein (CRP), Polyclonal Antibody Anti-Human Affinity Purified GCRP-80A

<b>Product Name</b>	C-Reactive Protein (CRP), Polyclonal Antibody Anti-Human Affinity Purified
<b>Catalog Number</b>	GCRP-80A
<b>Species/Host</b>	Goat
<b>Immunogen</b>	Highly purified human C-Reactive Protein (CRP)
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Clear solution
<b>Concentration</b>	1 - 2 mg/mL ( $A_{280}$ , $E_{1\text{mg/mL}}^{1\text{cm}}_{280\text{nm}} = 1.4$ )
<b>Storage</b>	2-8° C
<b>Purity</b>	Affinity Purified
<b>Source</b>	Polyclonal Antibody
<b>Shelf Life</b>	2 years