

Product Datasheet

C-Reactive Protein (CRP), Polyclonal Antibody Anti-Human Affinity Purified GCRP-80A

Product Name	C-Reactive Protein (CRP), Polyclonal Antibody Anti-Human Affinity Purified
Catalog Number	GCRP-80A
Species/Host	Goat
Immunogen	Highly purified human C-Reactive Protein (CRP)
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear solution
Concentration	1 - 2 mg/mL (A_{280} , $E_{1\text{mg/mL}}^{1\text{cm}}_{280\text{nm}} = 1.4$)
Storage	2-8° C
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years