

## Product Datasheet

### Retinol Binding Protein (RBP), Polyclonal Anti-Human Affinity Purified Host Goat GRBP-80A

<b>Product Name</b>	Retinol Binding Protein (RBP), Polyclonal Anti-Human Affinity Purified Host Goat
<b>Catalog Number</b>	GRBP-80A
<b>Species/Host</b>	Goat
<b>Immunogen</b>	High purity, native human RBP
<b>Alternative Names</b>	anti-RBP, anti-Retinol Binding Protein
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	2-8° C
<b>Isotype</b>	IgG
<b>Purity</b>	Affinity Purified
<b>Source</b>	Polyclonal Antibody
<b>Shelf Life</b>	2 years