

## Product Datasheet

### Retinol Binding Protein (RBP), Polyclonal Anti-Human Affinity Purified Host Goat GRBP-80A

|                          |  |
|--------------------------|--|
| <b>Product Name</b>      | Retinol Binding Protein (RBP), Polyclonal Anti-Human Affinity Purified Host Goat |
| <b>Catalog Number</b>    | GRBP-80A   |
| <b>Species/Host</b>      | Goat   |
| <b>Immunogen</b>         | High purity, native human RBP  |
| <b>Alternative Names</b> | anti-RBP, anti-Retinol Binding Protein   |
| <b>Brand</b>             | Lee Biosolutions   |
| <b>Form/Appearance</b>   | Liquid   |
| <b>Storage</b>           | 2-8° C   |
| <b>Isotype</b>           | IgG  |
| <b>Purity</b>            | Affinity Purified  |
| <b>Source</b>            | Polyclonal Antibody  |
| <b>Shelf Life</b>        | 2 years  |