

## Product Datasheet

### Anti-Human CD264 Monoclonal Antibody, Azide Free Clone B-L34 853.690

**Product Name** Anti-Human CD264 Monoclonal Antibody, Azide Free Clone B-L34

**Catalog Number** 853.690

**Description** Mouse monoclonal

**Species/Host** Mouse

**Conjugation** Azide Free

**Tested Applications** ELISA

**Immunogen** Recombinant human TRAIL R4/Fc chimera

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

**Isotype** Mouse IgG1 Kappa light chain

**Clonality** Monoclonal

**Specificity** Recognises TRAIL R4 (Dcr2), a 44.2 kDa protein

**Hybridoma** Myeloma X63/AG.8653 x Balb/c spleen cells