

Product Datasheet

Anti-Human L1CAM Monoclonal Antibody, Azide Free Clone B-R48 857.370

Product Name Anti-Human L1CAM Monoclonal Antibody, Azide Free Clone B-R48

Catalog Number 857.370

Species/Host Mouse

Conjugation Azide Free

Tested Applications ELISA

Immunogen Recombinant human L1CAM/Fc

Alternative Names CD171 / NCAM-L1

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Note Clone: B-R48

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognizes natural and recombinant human L1CAM

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells