

Product Datasheet

Anti-Human LRRC15 Monoclonal Antibody, Azide Free Clone B-G53 857.400

Product Name	Anti-Human LRRC15 Monoclonal Antibody, Azide Free Clone B-G53
Catalog Number	857.400
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	Recombinant human LRRC15
Alternative Names	Leucine rich repeat containing 15 / hLib / LIB protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE
Note	Clone: B-G53
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes natural and recombinant human LRRC15
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells