

Product Datasheet

Anti-Human LRRC15 Monoclonal Antibody, Unconjugated Clone B-G53 857.403.020

Product Name	Anti-Human LRRC15 Monoclonal Antibody, Unconjugated Clone B-G53
Catalog Number	857.403.020
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human LRRC15
Alternative Names	Leucine rich repeat containing 15 / hLib / LIB protein
Brand	Diaclone
Form/Appearance	Phosphate buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE
Note	Working Dilution: Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes natural and recombinant human LRRC15
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells