

Product Datasheet

Recombinant Human Leucine-Rich Repeat-Containing protein 15 LRRC15 (CHO) 705-H28-100

Product Name	Recombinant Human Leucine-Rich Repeat-Containing protein 15 LRRC15 (CHO)
Catalog Number	705-H28-100
Species/Host	CHO Cells
Alternative Names	LRRC15, LIB Protein, hLib
Brand	Diaclone
Form/Appearance	Liquid ; Phosphate-buffered saline pH7.4. Sterile-filtered through 0.22 µm.
Storage	Store at minimum -20°C. Avoid repeated freeze-thaw cycles.
Purity	>95% by SDS PAGE gel
Sequence	Met1-Gly538 (Extra Cellular Domain)
Tag	Human IgG1 Fc (Pro100-Lys330)