

## Product Datasheet

### Recombinant Human B Cell Maturation Protein (BCMA) 715-H29-100

<b>Product Name</b>	Recombinant Human B Cell Maturation Protein (BCMA)
<b>Catalog Number</b>	715-H29-100
<b>Species/Host</b>	HEK cells
<b>Alternative Names</b>	CD269, TNFRSF17, BCMA tumor necrosis factor receptor superfamily member 17
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Liquid. Phosphate-buffered saline pH7.4. Sterile-filtered through 0.22 µm
<b>Storage</b>	Store at minimum -20°C. Avoid repeated freeze-thaw cycles.
<b>Purity</b>	>90% by SDS PAGE gel
<b>Sequence</b>	ECD Leu2-Ala54
<b>Tag</b>	Silent Human IgG1 Fc (Pro100-Lys330)