

## Product Datasheet

### Anti-Human BCMA Monoclonal Antibody, Azide Free Clone B-B54 857.410

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Anti-Human BCMA Monoclonal Antibody, Azide Free Clone B-B54                                |
| <b>Catalog Number</b>      | 857.410  |
| <b>Species/Host</b>        | Mouse  |
| <b>Conjugation</b>         | Azide free   |
| <b>Tested Applications</b> | FC   |
| <b>Immunogen</b>           | Recombinant human BCMA ECD with Fc Human IgG1 domain                                       |
| <b>Alternative Names</b>   | CD269 / TNFRSF17 / B cell maturation antigen   |
| <b>Brand</b>               | Diaclone   |
| <b>Form/Appearance</b>     | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free |
| <b>Storage</b>             | Stable at +2-8°C for 12 months. For longer storage freeze aliquots.                        |
| <b>Note</b>                | Clone: B-B54   |
| <b>Isotype</b>             | Mouse IgG1 Kappa light chain   |
| <b>Source</b>              | Recombinant antibody from a phage display library  |
| <b>Specificity</b>         | Recognizes both natural and recombinant human BCMA   |