

## Product Datasheet

### Anti-Human BCMA Monoclonal Antibody, Unconjugated Clone B-B54 857.413.020

<b>Product Name</b>	Anti-Human BCMA Monoclonal Antibody, Unconjugated Clone B-B54
<b>Catalog Number</b>	857.413.020
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Recombinant human BCMA ECD with Fc Human IgG1 domain
<b>Alternative Names</b>	CD269 / TNFRSF17 / B cell maturation antigen
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate buffered saline with 1% BSA and 0.09% sodium azide
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
<b>Note</b>	Clone: B-B54
<b>Isotype</b>	Mouse IgG1 Kappa light chain
<b>Source</b>	Recombinant antibody from a phage display library
<b>Specificity</b>	Recognizes both natural and recombinant human BCMA