

Product Datasheet

Anti-Human CX3CL1 Monoclonal Antibody, Unconjugated Clone B-F51 857.423.020

Product Name	Anti-Human CX3CL1 Monoclonal Antibody, Unconjugated Clone B-F51
Catalog Number	857.423.020
Species/Host	Mouse
Conjugation	Purified
Tested Applications	ELISA, FC
Immunogen	Recombinant human CX3CL1
Alternative Names	Fractalkine / C-X3-C motif chemokine ligand 1
Brand	Diaclone
Form/Appearance	Phosphate buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE
Note	Clone: B-F51
Isotype	Mouse IgG2a Kappa light chain
Specificity	Recognize natural and recombinant human CX3CL1
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells