

## Product Datasheet

### Anti-Human CX3CL1 Monoclonal Antibody, Azide Free Clone B-K49 857.430

**Product Name** Anti-Human CX3CL1 Monoclonal Antibody, Azide Free Clone B-K49

**Catalog Number** 857.430

**Species/Host** Mouse

**Conjugation** Azide free

**Tested Applications** ELISA

**Immunogen** Recombinant Human CX3CL1 Full length

**Alternative Names** Fractalkine / C-X3-C motif chemokine ligand 1

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

**Note** Clone: B-K49

**Isotype** Mouse IgG1 Kappa light chain

**Specificity** Recognize natural and recombinant human CX3CL1

**Hybridoma** Myeloma X63/AG.8653 x Balb/c node cells