

Product Datasheet

Anti-Human CD48 Monoclonal Antibody, Azide Free Clone B-P59 857.440

Product Name Anti-Human CD48 Monoclonal Antibody, Azide Free Clone B-P59

Catalog Number 857.440

Species/Host Mouse

Conjugation Azide free

Tested Applications FC

Immunogen Recombinant human CD48

Alternative Names BLAST-1 / SLAMF2

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Note Clone: B-P59

Isotype Mouse IgG2a Kappa light chain

Specificity Recognizes natural and recombinant human CD48

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells