

Product Datasheet

Anti-Human CD48 Monoclonal Antibody, Unconjugated Clone B-P59 857.443.020

Product Name	Anti-Human CD48 Monoclonal Antibody, Unconjugated Clone B-P59
Catalog Number	857.443.020
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human CD48
Alternative Names	BLAST-1 / SLAMF2
Brand	Diaclone
Form/Appearance	Phosphate buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Clone: B-P59
Isotype	Mouse IgG2a Kappa light chain
Specificity	Recognizes natural and recombinant human CD48
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells