

## Product Datasheet

### Progesterone-3-HRP Antigen 170063

<b>Product Name</b>	Progesterone-3-HRP Antigen
<b>Catalog Number</b>	170063
<b>Description</b>	Synthetic Progesterone, HRP conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid; Antigen-HRP conjugate
<b>Storage</b>	Storage at 2-8 °C. Do not freeze.
<b>Note</b>	This product was previously available with a Catalog code P91-92-11H. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Tris base/FBS (1:1), pH 7.0, 0.0016% gentamicin sulfate, 0.005% thimerosal
<b>Matched Pairs</b>	100248, 100249, 100250, 100252, 100907, HM524, HM525
<b>Shelf Life</b>	Unspecified