

Product Datasheet

Benzoylecgonine (BE)-BSA Antigen 170048

Product Name	Benzoylecgonine (BE)-BSA Antigen
Catalog Number	170048
Description	Synthetic Benzoylecgonine, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short term storage (\leq 2 weeks): 2-8 °C. Long term:-20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P81-99-11A. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	100890 (detection)
Shelf Life	60 months