

## Product Datasheet

### Cotinine-4-BSA Antigen 170047

<b>Product Name</b>	Cotinine-4-BSA Antigen
<b>Catalog Number</b>	170047
<b>Description</b>	Synthetic Cotinine, BSA conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short term storage ( $\leq$ 2 weeks): 2-8 °C. Long term:-20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	This product was previously available with catalog code P82-99-04A. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Distilled water
<b>Matched Pairs</b>	PP194IF or PP198IF (detection)
<b>Shelf Life</b>	60 months