

## Product Datasheet

### BioCon II - Defibrinated Plasma H2000

<b>Product Name</b>	BioCon II - Defibrinated Plasma
<b>Catalog Number</b>	H2000
<b>Description</b>	BioCon II - Defibrinated Plasma
<b>Tested Applications</b>	QC Materials, Base Matrix
<b>Alternative Names</b>	Plasma Derived Serum, Filtered Plasma, Processed Plasma
<b>Brand</b>	Bioresource Technology
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Note</b>	This product is final formulated to the following parameters: Total Protein: >5.0 g/dL Sodium: 80-150 mmol/L Potassium: 2.0-6.0 mmol/L Chloride: 65-120 mmol/L Cholesterol: <200 mg/dL Triglycerides: <200 mg/dL pH: 7.0-7.6 Bioburden: <100 cfu/mL
<b>Purity</b>	0.2 micron filtered plasma
<b>Source</b>	Pooled Human Donors