

## Product Datasheet

### Cortisol-3-HRP Antigen 170053

<b>Product Name</b>	Cortisol-3-HRP Antigen
<b>Catalog Number</b>	170053
<b>Description</b>	Synthetic Cortisol
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid; Antigen-HRP conjugate
<b>Storage</b>	Storage at 2-8 °C. Do not freeze.
<b>Note</b>	This product was previously available with catalog code P91-92-91H. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Tris base/FBS (1:1), pH 7.0, 0.0016% gentamicin sulfate, 0.005% thimerosal
<b>Matched Pairs</b>	<b>Anti-Cortisol 12901 SPTN-5 100904,</b>
<b>Shelf Life</b>	Unspecified