

Product Datasheet

Anti-h NfL 12601 SPTN-5 100984

Product Name	Anti-h NfL 12601 SPTN-5
Catalog Number	100984
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	CLIA, FIA
Alternative Names	Neurofilament light chain
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Associated Products	NfL antigen LA666
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.6-7.3
Cross Reactivity	Does not recognize NfM, NfH, NSE or S100B. Recognizes bovine NfL.
Specificity	Antibody recognizes human neurofilament light
Matched Pairs	NfL 12603 and 12604
Shelf Life	Unspecified