

Product Datasheet

Anti-h NfL 12604 SPTN-5 100986

Product Name	Anti-h NfL 12604 SPTN-5
Catalog Number	100986
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	CLIA, FIA
Alternative Names	Neurofilament light chain
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 95 %
Associated Products	NfL antigen LA666
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.4-8.4
Cross Reactivity	Does not recognize NfM, NfH, NSE or S100B. Recognizes bovine NfL.
Specificity	Antibody recognizes human neurofilament light
Matched Pairs	NfL 12601
Shelf Life	Unspecified