

## Product Datasheet

### Immunoglobulin M (IgM) 340-39

<b>Product Name</b>	Immunoglobulin M (IgM)
<b>Catalog Number</b>	340-39
<b>Alternative Names</b>	IgM Concentrate, Partially Purified IgM
<b>Brand</b>	Lee Biosolutions
<b>Concentration</b>	≥ 15 mg/mL (Dimension® Clinical Chemistry System)
<b>Storage</b>	-20°C
<b>Purity</b>	Semi-pure
<b>Source</b>	Human Plasma (defibrinated) from Multiple Donors
<b>Buffer</b>	20 mM sodium phosphate, 0.5 M NaCl with 0.05% Proclin 300, pH 7.2,
<b>Shelf Life</b>	1 year