

## Product Datasheet

### Anti-Human CD117 Monoclonal Antibody, Azide Free Clone B-N20 853.730

**Product Name** Anti-Human CD117 Monoclonal Antibody, Azide Free Clone B-N20

**Catalog Number** 853.730

**Description** Mouse monoclonal

**Species/Host** Mouse

**Conjugation** Azide Free

**Tested Applications** ELISA

**Immunogen** Recombinant c-kit

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

**Isotype** Mouse IgG1 Kappa light chain

**Clonality** Monoclonal

**Specificity** Recognises the Stem Cell Factor Receptor, SCF R, c-kit, a 145 kDa protein

**Hybridoma** Myeloma X63/AG.8653 x Balb/c spleen cells