

Product Datasheet

Whole Blood, Pooled Human Donors 991-50-P

| | |
|------------------------|--|
| Product Name | Whole Blood, Pooled Human Donors |
| Catalog Number | 991-50-P |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid (Red) |
| Storage | 2-8°C |
| Note | Negative/non-reactive for HIV 1/2, HBsAg, RPR and HCV by FDA approved methods. Inquire for custom anticoagulants and donor pooling. |
| Source | Pooled Human Donors (≥ 3 donors per pool) |
| Shelf Life | Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time. |
| Preservatives | K3EDTA as an anticoagulant |