

Product Datasheet

Anti-h Gastrin-17 11701 SPTN-5 100788

Product Name	Anti-h Gastrin-17 11701 SPTN-5
Catalog Number	100788
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N-terminal part of gastrin-17
Purity	≥ 95 %
Associated Products	Gastrin 17 (G-17) Antigen LA451L
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.4-7.2
Cross Reactivity	Does not recognize gastrin-14 or gastrin-34
Specificity	Antibody recognizes human gastrin-17
Matched Pairs	Anti-h Gastrin-17 17005 and 17007
Shelf Life	12 months