

Product Datasheet

Anti-h Gastrin-17 11707 SPTN-5 100794

Product Name	Anti-h Gastrin-17 11707 SPTN-5
Catalog Number	100794
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG2a
Clonality	Monoclonal
Epitope	C-terminal part of gastrin-17
Purity	≥ 95 %
Associated Products	Gastrin 17 (G-17) Antigen LA451L
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.3-8.7
Cross Reactivity	Recognizes gastrin-14 and gastrin-34
Specificity	Antibody recognizes human gastrin-17
Matched Pairs	Anti-h Gastrin-17 17001 and 17002
Shelf Life	Unspecified