

Product Datasheet

Albumin (BSA, ALB) - Sulfhydryl Blocked 100-35

Product Name Albumin (BSA, ALB) - Sulfhydryl Blocked

Catalog Number 100-35

Alternative Names Serum ALB, HSA, ALB

Brand Lee Biosolutions

Form/Appearance Lyophilized

Concentration $\geq 98\%$ (Elemental analysis, dry basis)

Storage 2-8 °C

Purity $\geq 98\%$ (Electrophoresis)

Source Bovine Plasma

Molecular Weight 66,500