

## Product Datasheet

### Azurocidin 125-75

<b>Product Name</b>	Azurocidin
<b>Catalog Number</b>	125-75
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Lyophilized, essentially salt-free
<b>Concentration</b>	Reported in mg/mL ( $E^{0.1\%}_{280nm} = 0.51$ )
<b>Storage</b>	-20°C
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Human Neutrophils
<b>Molecular Weight</b>	37,000