

## Product Datasheet

### Beta-2 Microglobulin (B2M) 126-16

<b>Product Name</b>	Beta-2 Microglobulin (B2M)
<b>Catalog Number</b>	126-16
<b>Alternative Names</b>	B2M, BMG
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 - 4.0 mg B2M /mL (OD at 280nm)
<b>Storage</b>	2-8 °C
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Human Urine
<b>Buffer</b>	40 mM phosphate buffer, pH 7.4 ± 0.2 with 1 mM EDTA, 150 mM NaCl, and 0.1% sodium azide
<b>Molecular Weight</b>	~12,000