

Product Datasheet

Beta-2 Microglobulin (B2M) 126-16

Product Name	Beta-2 Microglobulin (B2M)
Catalog Number	126-16
Alternative Names	B2M, BMG
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	1.0 - 4.0 mg B2M /mL (OD at 280nm)
Storage	2-8 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Urine
Buffer	40 mM phosphate buffer, pH 7.4 ± 0.2 with 1 mM EDTA, 150 mM NaCl, and 0.1% sodium azide
Molecular Weight	~12,000