

Product Datasheet

Anti-Human CD263 Monoclonal Antibody, Azide Free Clone B-N29 853.740

Product Name	Anti-Human CD263 Monoclonal Antibody, Azide Free Clone B-N29
Catalog Number	853.740
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Recombinant human TRAIL R3/Fc chimera
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises Dcr1 (CD263), a 75.90kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells